

## Garp MRI Trolley SS

Type: Garp 00 With folding side rails.



**Intent of use:** For bringing a patient to the MRI scanner and slide the patient with over on MRI table.

### Features

- MRI Trolley 7T compatible
- Fiber tabletop 220 x 65cm made of antimagnetic material and painted off white
- X-Ray transparent tabletop.
- Handle in front and back.
- Hydraulic adjusted working height from 500 to 850mm
- Trendelenburg +/- 13,5 dg
- Adjustable backrest
- Side rails for accessories
- Folding side rails.
- Protective wheels in each corner
- Removable accessories tray
- Ø150mm wheels with foot brake, total and direction lock, all wheels can be controlled from one side.
- IV pole can be placed in each corner
- Two year of warranty
- Manufactured according to CE directive and CE-approval.
- Manufactured in Denmark