

LiGo MRI 3T Trolley EE



Type: G300515 0200



Inteded use:

MRI scanner: For bringing a patient to the MRI scanner to be scanned whilst on the transfer board

Features

- MRI compatible 3T
- MRI fiber transfer board with manual table locks
- Hydraulic powered - fully portable transfer board
- Hydraulic moved working height from 60 cm to 100cm
- Large wheels with central total and direction lock.
- Load capacity up to 250kg.

Technical data:

- Table frame is made of antimagnetic stainless steel and painted
- Distance between columns 74 cm
- MRI transfer board 200 x 56cm
- Special mattress 200 x 56 x 3cm
- Table frame, with side rails for mounting of side guards
- End and side rails for mounting of equipment.
- 150mm wheels with central triple-function foot brake, total lock, free and direction lock.
- Rubber brake pedals to avoid damage of ankles and legs
- Weight of table: 75kg
- Max. Weight of patient: 250kg.
- One year of warranty
- Manufactured according to CE directive with CE-approval
- Manufactured in Denmark

Product no. G300515 0200

OPTIONS

- Foldable side guards
- Drip stand
- Side handle
- Patient safety band
- Mattress 4 cm
- Arm support
- Patient side support
- Docking arm and stops