



**Knight
Imaging**

Cardiff Chair

A purpose built mammography chair developed with patient comfort, accessibility and operator use in mind.



- Fully adjustable chair for breast imaging, mammography and for stereotactic procedures.
- Designed to assist with imaging less mobile patients.
- Battery powered, motorised height adjustment.
- Manual quick release to assist patient recovery.

FEATURES

- Fully adjustable rear side sections and arm rests.
- Adjustable angle backrest from -5° to $+90^{\circ}$
- Manually adjustable foot section with height variable foot plate.
- Interchangeable head rest.
- Height Range: 60 cm - 100 cm
- Seat Width: 65 cm
- Seat Depth: 50 cm
- 10 cm thick comfort foam.
- Maximum Weight Capacity: 140 kg
- Anti-bacterial upholstery and paintwork.
- Castors can be set to steer front or rear, free movement, or locked in lateral pairs.



Paradigm House, 3 Melbourne Avenue, Broomhill, Sheffield, England, S10 2QJ, U.K.
Tel: +44 (0)114 267 0482, Fax: +44 (0)114 266 3944
E-mail: info@knightimaging.co.uk, Web: www.knightimaging.co.uk

A  LabLogic Company