



# Mammography Chair

[www.knightimaging.co.uk](http://www.knightimaging.co.uk)



**Knight  
Imaging**

EXPERIENCE & EXPERTISE

# Mammography Chair

EXPERIENCE & EXPERTISE

The VELA Mammography Chair is designed specifically for positioning during mammography examinations. Crafted to provide optimal working conditions for the mammographer this chair is a significant contributing factor to the production of higher quality images.

It provides good posture for staff, without unnecessary strain on the back, shoulders, and neck, since the work area can be adjusted to the mammographer regardless of their height.

Suitable for both CC and MLO views, as the patient can be placed and locked very close to the machine, and armrests can be folded down.

Optimal positioning of the patient on a secure and locked base, which, through the electric lift system, can easily be adjusted to the equipment and staff.

Greater image precision, as the patient can be secured in a forward-facing position, releasing a larger accessible breast area.



Electric seat height adjustment via a wired remote ensures the correct working height but also the ability to position the patient precisely at the X-ray machine.

Foot-operated brake and adjustable handlebar prevent inappropriate working postures.

**PATIENT WEIGHT**  
Up to 160 kgs  
(25 Stone)

Foldable armrests allow closer access to the machine during MLO imaging.

Easy to clean and disinfect

Smooth-running wheels make the chair easy to manoeuvre and position accurately at the machine.

Central brake ensures that the chair can be secured in a selected position until the examination is done and the mammographer is ready for repositioning.

Specifications	
Seat (W x D)	47 x 47 cm
Backrest (W x H)	40 x 40 cm
Seat Height	51 - 81 cm
Armrest Pad (W x L)	8 x 26 cm
Armrest Height (Measured from seat)	14 - 27 cm
Frame Size	55 x 55 cm

# Service and Support

Knight Imaging offers service and support on our range of couches and chairs. We carry out one-off service calls, yearly preventative maintenance visits as well as offering support contracts.

## Planned Preventative Maintenance

Our fully trained engineers service your couches and chairs ensuring that any potential issues are spotted and brought to your attention before they become a problem.

We typically carry out these types of services across entire departments ensuring the best possible lifetime of your couches and chairs. We will always arrange visits in advance to ensure a swift turnaround and to minimise disruption to your ward or department.

A preventative maintenance visit can be arranged whenever you think you need to have your equipment inspected and checked for faults or potential future problems. Just get in touch to discuss your equipment and make a booking.

## Breakdown Response

Our experienced engineers are available to diagnose and repair in times of emergency.



**Knight Imaging Limited**  
Innovation House, 6 Europa View  
Sheffield, S9 1XH, UK  
**E-mail:** [sales@knightimaging.co.uk](mailto:sales@knightimaging.co.uk)  
**Tel:** +44 (0)114 267 0482  
**[www.knightimaging.co.uk](http://www.knightimaging.co.uk)**

