



Made in
England

Autotilt Sigma

Specialist Critical Care Chairs

www.knightimaging.co.uk



**Knight
Imaging**

EXPERIENCE & EXPERTISE

Specialist chairs for critical and bariatric care

The range of battery powered Sigmas, which are designed to meet patient needs in Critical Care units, are ideal for enabling staff to move a patient out of bed at the earliest opportunity to aid recovery.

From the operator's perspective the powered controls make for easier movement and optimum positioning. In addition the Quick Release function gives staff the peace of mind that they can bring the patient into a horizontal position in under 3 seconds.

Tilt in space

The battery motors operate the variable height, chair to trolley mode and Trendelenburg positioning to +/- 10°. The leg rests are independently adjustable for angulation and length. Arm rests are height adjustable.

Quick release

Manual quick release facility for horizontal positioning in seconds.



Choice of models

Choice of models available for patients small or large, slight or obese.

A range of accessories

Drip Stand

Double hook, 2 section adjustable height. Can be fitted to either side of the Sigma.



Oxygen Bottle Holder

Fits white CD2 O₂ cylinder which can be fitted vertically to either side of the Sigma.



Lateral Supports

Adjustable inboard reach with matching anti-bacterial upholstery to the Sigma. 2 positions: Upper and lower or backrest.



Castors can be set to steer front or rear, free movement, or locked in lateral pairs

Autotilt Sigma 315 – Basic Specifications

| | |
|---|-------------------|
| Maximum Patient Weight Capacity | 135 kg / 21 stone |
| Chair Width / Depth (Inside) | 54 cm / 44 cm |
| Chair Mode Height Range | 65 - 86 cm |
| Trolley Mode Height Range | 84 - 106 cm |
| Please refer to the Technical Specification Sheet for further information | |

Autotilt Sigma 250 – Basic Specifications

| | |
|---|-------------------|
| Maximum Patient Weight Capacity | 250 kg / 39 stone |
| Chair Width / Depth (Inside) | 74 cm / 44 cm |
| Chair Mode Height Range | 64 - 84 cm |
| Trolley Mode Height Range | 82 - 103 cm |
| Please refer to the Technical Specification Sheet for further information | |

Autotilt Sigma 300 – Basic Specifications

| | |
|---|-------------------|
| Maximum Patient Weight Capacity | 300 kg / 47 stone |
| Chair Width / Depth (Inside) | 74 cm / 54 cm |
| Chair Mode Height Range | 66 - 86 cm |
| Trolley Mode Height Range | 84 - 106 cm |
| Please refer to the Technical Specification Sheet for further information | |

Service and Support

Knight Imaging offers service and support on our range of couches and chairs. We carry out one-off service calls, yearly preventative maintenance visits as well as offering support contracts.

Planned Preventative Maintenance

Our fully trained engineers service your couches and chairs ensuring that any potential issues are spotted and brought to your attention before they become a problem.

We typically carry out these types of services across entire departments ensuring the best possible lifetime of your couches and chairs. We will always arrange visits in advance to ensure a swift turnaround and to minimise disruption to your ward or department.

A preventative maintenance visit can be arranged whenever you think you need to have your equipment inspected and checked for faults or potential future problems. Just get in touch to discuss your equipment and make a booking.

Breakdown Response

Fully prepared, with a team of experienced engineers and an availability of parts for all the couches and chairs in the range, Knight Imaging is ideally placed to diagnose and repair in times of emergency.

Knight Imaging Limited

Paradigm House, 3 Melbourne Avenue
Broomhill, Sheffield, S10 2QJ, UK

E-mail: sales@knightimaging.co.uk

Tel: +44 (0)114 267 0482

Fax: +44 (0)114 266 3944

www.knightimaging.co.uk

