



Made in  
England

# Cardiff Chair

Specialist Mammography Chair

[www.knightimaging.co.uk](http://www.knightimaging.co.uk)



**Knight  
Imaging**

EXPERIENCE & EXPERTISE

## A purpose built mammography chair developed with patient comfort, accessibility and operator use in mind

A fully adjustable chair for breast imaging, mammography and for stereotactic procedures.

### Battery powered

Designed to assist with imaging less mobile patients, the Cardiff Chair is battery powered with motorised height adjustment.

### Quick release

Manual quick release in case of emergency, as well as to assist patient recovery.



#### Basic Specifications

Maximum Patient Weight Capacity	140 kg / 22 stone
Height Range	60 - 100 cm
Seat Width	65 cm
Seat Depth	

Please refer to the Technical Specification Sheet for further information

# Service and Support

Knight Imaging offers service and support on our range of couches and chairs. We carry out one-off service calls, yearly preventative maintenance visits as well as offering support contracts.

## Planned Preventative Maintenance

Our fully trained engineers service your couches and chairs ensuring that any potential issues are spotted and brought to your attention before they become a problem.

We typically carry out these types of services across entire departments ensuring the best possible lifetime of your couches and chairs. We will always arrange visits in advance to ensure a swift turnaround and to minimise disruption to your ward or department.

A preventative maintenance visit can be arranged whenever you think you need to have your equipment inspected and checked for faults or potential future problems. Just get in touch to discuss your equipment and make a booking.

## Breakdown Response

Fully prepared, with a team of experienced engineers and an availability of parts for all the couches and chairs in the range, Knight Imaging is ideally placed to diagnose and repair in times of emergency.

### **Knight Imaging Limited**

Paradigm House, 3 Melbourne Avenue  
Broomhill, Sheffield, S10 2QJ, UK

**E-mail:** [sales@knightimaging.co.uk](mailto:sales@knightimaging.co.uk)

**Tel:** +44 (0)114 267 0482

**Fax:** +44 (0)114 266 3944

[www.knightimaging.co.uk](http://www.knightimaging.co.uk)

