



Cabinet Trolley - TR149

Designed for use with electrotherapy devices, this mobile steel cabinet provides both storage and security with 3 internal shelves for positioning the treatment units and accessories.

The trolley and doors are powder coated to provide a hard wearing and easy to clean finish.



- Fully welded construction making this a strong and durable trolley
- 75mm twin wheel easy glide castors
- Non marking castors

INFORMATION

Dimensions: L400mm x W400mm x H980mm

Top shelf: 400mm x 300mm

Internal shelves: 315mm x 210mm

Freight Dimensions:

450mm x 450mm x 900mm

Weight: 21kg