



Operators Chairs and Stools

Deluxe Operators Chair



Scanning Chair



High Operators Stool



Deluxe Saddle Seat



High Operators Chair



Standard Operators Chair

ORDER CODE	DESCRIPTION
M7A-NXX-01	Standard Operators Chair <ul style="list-style-type: none"> Height range 48 - 60 cm (19" - 24") with adjustable backrest. 5 star base with castors.
M7A-NXX-02	Standard Operators Chair <ul style="list-style-type: none"> Auto-step backrest. Seat tilt. Height range 48 - 60 cm (19" - 24"). 5 star base with castors.
M7A-NXX-03	High Operators Chair <ul style="list-style-type: none"> Height range 59 - 84 cm (23" - 33") with adjustable backrest. Adjustable height foot ring support. 5 star base with glides.
M7A-NXX-04	High Operators Chair <ul style="list-style-type: none"> Auto-step backrest. Seat tilt. Height range 59 - 84 cm (23" - 33"). Adjustable height foot ring support. 5 star base with glides.
M7A-NXX-05	Deluxe Operators Chair <ul style="list-style-type: none"> Height range 48 - 60 cm (19" - 24") with adjustable high level backrest. Adjustable height armrest supports (5 settings). Seat tilt. 5 star base with castors.
M7A-NXX-06	Deluxe Operators Chair <ul style="list-style-type: none"> Height range 48 - 60 cm (19" - 24") with auto-step backrest. Adjustable height armrest supports (5 settings). Seat tilt. 5 star base with castors.
M7A-NXX-07	Standard Operators Stool <ul style="list-style-type: none"> Height range 48 - 60 cm (19" - 24") 5 star base with castors.
M7A-NXX-08	High Operators Stool <ul style="list-style-type: none"> Height range 59 - 84 cm (23" - 33"). Adjustable foot ring support. 5 star base with glides.

ORDER CODE	DESCRIPTION
M7A-NXX-09	Mobile Therapy Stool <ul style="list-style-type: none"> Height range 47 - 60 cm. Standard model with castors.
M7A-NXX-10	Therapy Saddle Chair <ul style="list-style-type: none"> Standard model with seat tilt and castors. Height range 54 - 71 cm.
M7A-NXX-11	Deluxe Saddle Seat <ul style="list-style-type: none"> Tilting seat facility. Height range 54 - 71 cm (63 - 84 cm available on request). 5 star base with castors.
M7A-NXX-12	Therapy Saddle Chair <ul style="list-style-type: none"> Deluxe model with seat tilt and castors. Height range 63 - 84 cm.
M7A-NXX-13	Economy Saddle Seat <ul style="list-style-type: none"> Height range 48 - 60 cm (19" - 24") with adjustable high level backrest. 5 star base with castors.
M7A-NXX-14	Operators Work Stool <ul style="list-style-type: none"> PU model with foot ring. Height 47 - 60 cm.
M7A-NXX-15	Operators Work Chair <ul style="list-style-type: none"> PU model with backrest. Height 47 - 60 cm.
M7A-NXX-16	Operators Work Chair <ul style="list-style-type: none"> High PU model with backrest and glides.
M7A-NXX-17	Scanning Chair <ul style="list-style-type: none"> Spring loaded backrest support and base with castors. Height 48 - 60 cm (19" - 24").



Standard Operators Stool

Specifications are subject to change without notice.



Knight Imaging Limited

Paradigm House, 3 Melbourne Avenue
Broomhill, Sheffield, S10 2QJ, UK

E-mail: sales@knightimaging.co.uk

Tel: +44 (0)114 267 0482

Fax: +44 (0)114 266 3944

www.knightimaging.co.uk

Version 1.0 January 2019